




Product specification

ESR rack, suitable for Microvette®, 10 measuring stations



Product information

Product number	90.1091
----------------	---------

Product descriptions

Name	ESR rack, suitable for Microvette®, 10 measuring stations
Product number	90.1091
Product description	ESR rack, suitable for Microvette®, 10 measuring stations

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Product characteristics

Number of positions	10
---------------------	----

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Piece(s) / pallet	160
Height of case	335 mm
Width of case	98 mm

Depth of case	263 mm
Volume of case	0.0086 cbm
Weight of case	5.1 kg
EAN inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

