




## Product specification

### ESR rack, suitable for Microvette®, 10 measuring stations



## Product information

Product number	90.1091
----------------	---------

## Product descriptions

Name	ESR rack, suitable for Microvette®, 10 measuring stations
Product number	90.1091
Product description	ESR rack, suitable for Microvette®, 10 measuring stations

## Product characteristics

Number of positions	10
---------------------	----

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

## Packaging

Piece(s) / pallet	160
Height of inner box	20 mm
Width of inner box	260 mm
Depth of inner box	305 mm
Height of case	340 mm

Width of case	400 mm
Depth of case	585 mm
Volume of case	0.0796 cbm
Weight of case	5.1 kg
EAN inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT