



## Product specification

### Microvette® 300 Lithium heparin LH, 300 µl, cap green, flat base



## Product information

Product number	20.1309.100
----------------	-------------

## Product descriptions

Name	Microvette® 300 Lithium heparin LH, 300 µl, cap green, flat base
------	--

Product number	20.1309.100
----------------	-------------

Product description	Microvette® 300 Lithium heparin LH, capillary blood collection, preparation: lithium heparin LH, nominal volume: 300 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap green, flat base, inner receptacle: conical, 100 piece(s)/case
---------------------	---

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

## Size

Diameter	10.8 mm
Length excluding cap	46.6 mm
Length including cap	47.6 mm

Sample volume	300 µl
---------------	--------

## Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	green
Shape of inner tube	conical
Cap	cap
Type of collection	capillary
Centrifugation recommended	2000 U/min
Centrifugation time	5 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	121000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	2.43 kg
EAN inner box	4038917127615
EAN case	4038917021296

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT