



Product specification

Blood gas capillaries, calcium-balanced lithium heparin, 175 µl



Product information

| | |
|----------------|------------|
| Product number | 19.930.175 |
|----------------|------------|

Product descriptions

| | |
|---------------------|---|
| Name | Blood gas capillaries, calcium-balanced lithium heparin, 175 µl |
| Product number | 19.930.175 |
| Product description | Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 175 µl, (LxØ): 100 x 2.3 mm, 200 piece(s)/container |

Material & colours

| | |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

Size

| | |
|-------------------|--------|
| Diameter | 2.3 mm |
| Length of product | 100 mm |
| Sample volume | 175 µl |

Product characteristics

| | |
|---------------------------|----------------|
| Preparation concentration | 80 IU/ml blood |
| Type of collection | capillary |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Latex-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------|---------------|
| Type of inner packaging | Case |
| Piece(s) / inner box | 200 |
| Piece(s) / case | 12800 |
| Piece(s) / pallet | 921600 |
| Height of case | 190 mm |
| Width of case | 303 mm |
| Depth of case | 333 mm |
| Volume of case | 0.0192 cbm |
| Weight of case | 5.504 kg |
| EAN inner box | 4038917234344 |
| EAN case | 4038917234337 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT