




Product specification

Minivette® POCT EDTA K3E, 200 µl, plunger violet, colour code ISO, 150 piece(s)/bag



Product information

Product number	17.2113.201
----------------	-------------

Product descriptions

Name	Minivette® POCT EDTA K3E, 200 µl, plunger violet, colour code ISO, 150 piece(s)/bag
Product number	17.2113.201
Product description	Minivette® POCT EDTA K3E, nominal volume: 200 µl, preparation: K3 EDTA, violet, colour code ISO, 150 piece(s)/bag, 150 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Sample volume	200 µl
---------------	--------

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration

Latex-free	yes
Batched	yes

Packaging

Piece(s) / inner box	150
Piece(s) / case	1500
Piece(s) / pallet	108000
Height of inner box	117 mm
Width of inner box	99 mm
Depth of inner box	140 mm
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm
Weight of case	2.3595 kg
EAN inner box	4038917362177
EAN case	4038917362160

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT