



## Product specification

**Minivette® POCT EDTA K3E, 100 µl, plunger violet, colour code ISO, 200 piece(s)/bag**



## Product information

|                |             |
|----------------|-------------|
| Product number | 17.2113.101 |
|----------------|-------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Minivette® POCT EDTA K3E, 100 µl, plunger violet, colour code ISO, 200 piece(s)/bag  |
| Product number      | 17.2113.101  |
| Product description | Minivette® POCT EDTA K3E, nominal volume: 100 µl, preparation: K3 EDTA, violet, colour code ISO, 200 piece(s)/bag, 200 piece(s)/case |

## Material & colours

|                  |                        |
|------------------|------------------------|
| Product material | Cellulose acetate (CA) |
| Color of product | transparent            |

## Size

|               |        |
|---------------|--------|
| Sample volume | 100 µl |
|---------------|--------|

## Product characteristics

|                    |           |
|--------------------|-----------|
| Type of collection | capillary |
|--------------------|-----------|

## Purity & certification

|                  |                       |
|------------------|-----------------------|
| Product category | diagnostic device, CE |
|------------------|-----------------------|

|              |     |
|--------------|-----|
| CE certified | CE  |
| Latex-free   | yes |
| Batched      | yes |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 2000          |
| Piece(s) / inner box   | 200           |
| Piece(s) / case        | 2000          |
| Piece(s) / pallet      | 144000        |
| Height of inner box    | 117 mm        |
| Width of inner box     | 99 mm         |
| Depth of inner box     | 140 mm        |
| Height of case         | 140 mm        |
| Width of case          | 293 mm        |
| Depth of case          | 513 mm        |
| Volume of case         | 0.021 cbm     |
| Weight of case         | 2.51 kg       |
| EAN inner box          | 4038917375597 |
| EAN case               | 4038917375580 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT