



Product specification

Microvette® CB 300 EDTA K2E, 300 µl, cap violet, push cap, flat base



Product information

Product number	16.444.100
----------------	------------

Product descriptions

Name	Microvette® CB 300 EDTA K2E, 300 µl, cap violet, push cap, flat base
Product number	16.444.100
Product description	Microvette® CB 300 EDTA K2E, Veterinary medicine, preparation: K2 EDTA, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap violet, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	violet

Size

Diameter	10.8 mm
Length excluding cap	43.7 mm
Length including cap	46.5 mm
Sample volume	300 µl

Product characteristics

Base shape	Flat base
Shape of inner tube	trumpet-shaped
Closure type	Push cap
Type of collection	capillary

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Height of case	180 mm
Width of case	288 mm
Depth of case	378 mm
Volume of case	0.0196 cbm
Weight of case	2.39 kg
EAN inner box	4038917268356
EAN case	4038917282772

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT