



Product specification

Haemo-Diff®, with smear edge, for producing blood smears



Product information

Product number	14.1218
----------------	---------

Product descriptions

Name	Haemo-Diff®, with smear edge, for producing blood smears
Product number	14.1218
Product description	Haemo-Diff®, with smear edge, for producing blood smears, 500 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	21 mm
-------------------	-------

Product characteristics

Design	with smear edge
Connection	S-Monovette

Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	216000
Height of inner box	129 mm
Width of inner box	129 mm
Depth of inner box	191 mm
Height of case	150 mm
Width of case	268 mm
Depth of case	398 mm
Volume of case	0.016 cbm
Weight of case	3.27 kg
EAN inner box	4038917085830
EAN case	4038917004466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT