



## Product specification

### Haemo-Diff®, with smear edge, for producing blood smears



## Product information

Product number	14.1218
----------------	---------

## Product descriptions

Name	Haemo-Diff®, with smear edge, for producing blood smears
------	--

Product number	14.1218
----------------	---------

Product description	Haemo-Diff®, with smear edge, for producing blood smears, 500 piece(s)/case
---------------------	---

## Material & colours

Product material	Polystyrene (PS)
------------------	------------------

Color of product	transparent
------------------	-------------

## Size

Length of product	21 mm
-------------------	-------

## Product characteristics

Design	with smear edge
--------	-----------------

Connection	S-Monovette
------------	-------------

## Purity & certification

Product category	no medical device   IVD
------------------	-------------------------

Batched

yes

## Packaging

Minimum order quantity	2000
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	216000
Height of inner box	129 mm
Width of inner box	129 mm
Depth of inner box	191 mm
Height of case	150 mm
Width of case	268 mm
Depth of case	398 mm
Volume of case	0.016 cbm
Weight of case	3.27 kg
EAN inner box	4038917085830
EAN case	4038917004466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT