



Product specification

Haemo-Diff®, without smear edge, for producing blood smears



Product information

Product number	14.1217
----------------	---------

Product descriptions

Name	Haemo-Diff®, without smear edge, for producing blood smears
Product number	14.1217
Product description	Haemo-Diff®, without smear edge, for producing blood smears, 500 piece(s)/case

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	26.5 mm
-------------------	---------

Product characteristics

Design	without smear edge
--------	--------------------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	648000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.81 kg
EAN inner box	4038917085823
EAN case	4038917004459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT