



# Product specification

## Tourniquet, latex-free



### Product information

Product number	95.1007.101
----------------	-------------

### Product descriptions

Name	Tourniquet, latex-free
Product number	95.1007.101
Product description	Tourniquet, latex-free, 250 piece(s)/bag

### Material & colours

Color of product	yellow
------------------	--------

### Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	yes

### Packaging

Piece(s) / inner box	250
Piece(s) / case	2000
Piece(s) / pallet	42000
Volume of case	0 cbm
Weight of case	7.01 kg
EAN inner box	4038917263214

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT