



# Product specification

Carrier tube, (LxØ): 60 x 11.5 mm, PP, brown



## Product information

Product number	55.1577.002
----------------	-------------

## Product descriptions

Name	Carrier tube, (LxØ): 60 x 11.5 mm, PP, brown
Product number	55.1577.002
Product description	Carrier tube, (LxØ): 60 x 11.5 mm, material: PP, brown, 250 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	brown

## Size

Diameter	11.5 mm
Length of product	60 mm

## Product characteristics

Base shape	Flat base
Light protection	yes

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	78000
Height of case	130 mm
Width of case	323 mm
Depth of case	428 mm
Volume of case	0.018 cbm
Weight of case	1.7 kg
EAN inner box	4038917131032
EAN case	4038917023597

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT