



## Product specification

### Plunger, translucent, for Seraplas® V filter



## Product information

Product number	53.1030
----------------	---------

## Product descriptions

Name	Plunger, translucent, for Seraplas® V filter
Product number	53.1030
Product description	Plunger, translucent, for Seraplas® V filter, 200 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	translucent

## Size

Length of product	89 mm
-------------------	-------

## Purity & certification

Product category	no medical device   IVD
Batched	no
Safety data sheet	no

## Packaging

Minimum order quantity	2000
Piece(s) / inner box	200

Piece(s) / case	2000
Piece(s) / pallet	96000
Height of case	215 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0315 cbm
Weight of case	3.12 kg
EAN inner box	4038917085076
EAN case	4038917003773

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT