



## Product specification

**Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugat**



## Product information

Product number	53.677
----------------	--------

## Product descriptions

Name	Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation
------	---

Product number	53.677
----------------	--------

Product description	Seraplas® valve filter, V11, white, for separation of serum/plasma from the blood cells after centrifugation, 500 piece(s)/bag
---------------------	--

## Material & colours

Color of product	white
------------------	-------

## Size

Length of product	92 mm
-------------------	-------

## Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Batched	no
---------	----

Safety data sheet	no
-------------------	----

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of inner box	75 mm
Width of inner box	300 mm
Depth of inner box	300 mm
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	8.5 kg
EAN inner box	4038917083317
EAN case	4038917002257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT