



Product specification

Quick-release cap 1.60–2.10 mm for blood gas capillaries

140 µl / 175 µl, blue



Product information

Product number	65.935.230
----------------	------------

Product descriptions

Name	Quick-release cap 1.60–2.10 mm for blood gas capillaries 140 µl / 175 µl, blue
Product number	65.935.230
Product description	Quick-release cap 1.60–2.10 mm for blood gas capillaries 140 µl / 175 µl, blue, blood gas analysis, 400 piece(s)/bag

Material & colours

Color of product	blue
------------------	------

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

Packaging

Piece(s) / inner box	400
Piece(s) / case	32000
Piece(s) / pallet	2048000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	8.128 kg
EAN inner box	4038917234580
EAN case	4038917234573

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT