



## Product specification

**Ventilator for Blood Gas Monovette®, transparent/orange,  
length: 24.5 mm**



## Product information

Product number	14.1148
----------------	---------

## Product descriptions

Name	Ventilator for Blood Gas Monovette®, transparent/orange, length: 24.5 mm
Product number	14.1148
Product description	Ventilator for Blood Gas Monovette®, transparent/orange, length: 24.5 mm, accessories for the Blood Gas Monovette®, 500 piece(s)/case

## Material & colours

Color of product	transparent/orange
------------------	--------------------

## Size

Length of product	24.5 mm
-------------------	---------

## Product characteristics

Example of use	Ventilation of the Blood Gas - Monovette®
----------------	---

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration

Latex-free	yes
Batched	yes

## Packaging

Type of inner packaging	Case
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	168000
Height of inner box	129 mm
Width of inner box	129 mm
Depth of inner box	191 mm
Height of case	150 mm
Width of case	218 mm
Depth of case	288 mm
Volume of case	0.0094 cbm
Weight of case	1.34 kg
EAN inner box	4038917085465
EAN case	4038917075312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT