


Product specification

Membrane adapter



Product information

| | |
|--------------|---------|
| Order number | 14.1112 |
|--------------|---------|

Product descriptions

| | |
|------|------------------|
| Name | Membrane adapter |
|------|------------------|

| | |
|----------------|---------|
| Product number | 14.1112 |
|----------------|---------|

| | |
|---------------------|---|
| Product description | Membrane adapter, connection: Luer, for S-Monovette® needle/Multify® needle, 1 piece(s)/blister, 100 piece(s)/case, sterile |
|---------------------|---|

Product characteristics

| | |
|------------|------|
| Connection | Luer |
|------------|------|

| | |
|----------------|--|
| Example of use | Adaptation of an S-Monovette® Needle or (Safety) Multify® Needle to a Luer system, e.g. Blood-Gas Monovette® |
|----------------|--|

| | |
|--------------|-------------------------------------|
| Suitable for | S-Monovette® needle/Multify® needle |
|--------------|-------------------------------------|

Size

| | |
|-------------------|-------|
| Length of product | 16 mm |
|-------------------|-------|

Material & colours

| | |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
|------------------|-----------------------------------|

| | |
|------------------|-------------|
| Color of product | translucent |
|------------------|-------------|

Purity & certification

| | |
|-----------------------------|-------------------------|
| Accordance with MDD or IVDD | In-vitro diagnostic, CE |
|-----------------------------|-------------------------|

| | |
|--------------|----|
| CE certified | CE |
|--------------|----|

| | |
|-----------------|---------|
| Purity standard | sterile |
|-----------------|---------|

| | |
|---------------|----------------------|
| Sterilisation | Electron irradiation |
|---------------|----------------------|

| | |
|--------------|-----|
| Pyrogen-free | yes |
|--------------|-----|

| | |
|----------------------|-----|
| Lot or serial number | yes |
|----------------------|-----|

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 100 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 80000 |
| Height of inner box | 71 mm |
| Width of inner box | 81 mm |
| Depth of inner box | 268 mm |
| Height of case | 287 mm |
| Width of case | 178 mm |
| Depth of case | 373 mm |
| Volume of case | 0.0191 cbm |
| Weight of case | 1.49 kg |
| EAN inner box | 4038917085359 |
| EAN case | 4038917004053 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

