




Product specification

Membrane adapter



Product information

Product number	14.1112
----------------	---------

Product descriptions

Name	Membrane adapter
Product number	14.1112
Product description	Membrane adapter, connection: Luer, for S-Monovette® needle/Multifly® needle, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

Material & colours

Color of product	translucent
------------------	-------------

Size

Length of product	16 mm
-------------------	-------

Product characteristics

Connection	Luer
Example of use	Adaptation of an S-Monovette® Needle or (Safety) Multifly® Needle to a Luer system, e.g. Blood-Gas Monovette®

Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Height of inner box	71 mm
Width of inner box	81 mm
Depth of inner box	268 mm
Height of case	178 mm
Width of case	287 mm
Depth of case	373 mm
Volume of case	0.0191 cbm
Weight of case	1.49 kg
EAN inner box	4038917085359
EAN case	4038917004053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT