




Product specification

**Safety-Multifly® needle, 21G x 3/4", green, tube length:
200 mm, 1 piece(s)/blister**



Product information

Product number	85.1638.935
----------------	-------------

Product descriptions

Name	Safety-Multifly® needle, 21G x 3/4", green, Tube length: 200 mm, 1 piece(s)/blister
Product number	85.1638.935
Product description	Safety-Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 200 mm, with needle holder assembled, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 25 piece(s)/case

Size

Diameter of tube	2.3 mm
Length of tube	200 mm

Product characteristics

Design	needle holder assembled
Safety product	yes
Connection	Vacuum tube

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	7200
Height of case	140 mm
Width of case	263 mm
Depth of case	533 mm
Volume of case	0.0196 cbm
Weight of case	1.38 kg
EAN inner box	4038917240345
EAN case	4038917240338

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT