




## Product specification

**Multifly® needle, 21G x 3/4", green, tube length: 240 mm,  
1 piece(s)/blister**



## Product information

Product number	85.1638.001
----------------	-------------

## Product descriptions

Name	Multifly® needle, 21G x 3/4", green, Tube length: 240 mm, 1 piece(s)/blister
Product number	85.1638.001
Product description	Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 240 mm, with Luer (f), DEHP-free, 1 piece(s)/blister, sterile, non- pyrogenic/endotoxin-free, 120 piece(s)/case

## Size

Diameter of tube	2.3 mm
Length of tube	240 mm

## Product characteristics

Design	Luer (f)
Safety product	no
Connection	Luer

## Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

## Packaging

Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	20640
Height of inner box	171 mm
Width of inner box	214 mm
Depth of inner box	224 mm
Height of case	240 mm
Width of case	358 mm
Depth of case	442 mm
Volume of case	0.038 cbm
Weight of case	2.9198 kg
EAN inner box	4038917423755
EAN case	4038917423465

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT