



Product specification

S-Monovette® safety needle, 22G x 1 1/2", black, 1 piece(s)/blister



Product information

Product number	85.1440.200
----------------	-------------

Product descriptions

Name	S-Monovette® safety needle, 22G x 1 1/2", black, 1 piece(s)/blister
------	---------------------------------------------------------------------

Product number	85.1440.200
----------------	-------------

Product description	S-Monovette® safety needle, needle (ØxL): 22G x 1 1/2", coding colour: black, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case
---------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------

Product characteristics

Safety product	yes
----------------	-----

Connection	Membrane screw cap
------------	--------------------

Purity & certification

Product category	Medical device (MD)
------------------	---------------------

CE certified	CE0124 - DEKRA
--------------	----------------

Purity standard	sterile, non-pyrogenic/endotoxin-free
-----------------	---------------------------------------

Sterilisation	Gamma irradiation
---------------	-------------------

Pyrogen-free	yes
--------------	-----

Batched

yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Height of inner box	71 mm
Width of inner box	111 mm
Depth of inner box	268 mm
Height of case	290 mm
Width of case	233 mm
Depth of case	373 mm
Volume of case	0.0252 cbm
Weight of case	2.47 kg
EAN inner box	4038917129190
EAN case	4038917361965

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT