



## Product specification

### Multi-adapter, Luer lock



## Product information

Product number	14.1205.100
----------------	-------------

## Product descriptions

Name	Multi-adapter, Luer lock
Product number	14.1205.100
Product description	Multi-adapter, Luer lock, connection: Luer lock, for S-Monovette®, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

## Material & colours

Color of product	translucent
------------------	-------------

## Size

Length of product	47 mm
-------------------	-------

## Product characteristics

Design	Luer lock
Connection	Luer lock

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation

Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Height of inner box	71 mm
Width of inner box	81 mm
Depth of inner box	268 mm
Height of case	178 mm
Width of case	287 mm
Depth of case	373 mm
Volume of case	0.0191 cbm
Weight of case	3.11 kg
EAN inner box	4038917126700
EAN case	4038917020817

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT