



Product specification

Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled



Product information

Product number	14.1209
----------------	---------

Product descriptions

Name	Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled
Product number	14.1209
Product description	Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 20 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	84 mm
Diameter of opening	38 mm

Product characteristics

Design	Membrane adapter assembled
--------	----------------------------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	20
Piece(s) / case	100
Piece(s) / pallet	3600
Height of inner box	170 mm
Width of inner box	110 mm
Depth of inner box	350 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	1.94 kg
EAN inner box	4038917354776
EAN case	4038917354769

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT