



Product specification

Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled



Product information

Product number	14.1207
----------------	---------

Product descriptions

Name	Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled
Product number	14.1207
Product description	Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 50 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	74 mm
Diameter of opening	19 mm

Product characteristics

Design	membrane screw cap assembled
--------	------------------------------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	9000
Height of inner box	170 mm
Width of inner box	110 mm
Depth of inner box	350 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	2.14 kg
EAN inner box	4038917364188
EAN case	4038917364171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT